

Fecal Float is a standardized solution of sodium nitrate with a specific gravity of 1.25 to 1.30.

INDICATIONS: For use in all flotation techniques of fecal analysis.

Iss. 09-12

NDC 57319-605-09

Fecal Float Ready To Use

Sodium Nitrate

**FOR ANIMAL USE ONLY
KEEP OUT OF REACH OF CHILDREN**

ACTIVE INGREDIENT: Sodium Nitrate

INACTIVE INGREDIENT: Water

STORAGE: Store at controlled room temperature between 15°-30°C (59°-86°F). Keep container tightly closed when not in use.

Protect from freezing.

SHAKE WELL BEFORE USING.



Lot No.

Exp. Date

Trademarks are property of
Clipper Distributing Company, LLC



**Net Contents:
One Gallon (3.785 L)**

**PHOENIX™**
Manufactured for:
Clipper Distributing Company, LLC
St. Joseph, MO 64507

Manufactured by:
First Priority, Inc.
Elgin, IL U.S.A.